

The logo for SmitCare, featuring the word "SmitCare" in a sans-serif font. "Smit" is in black and "Care" is in red. The text is underlined with a red line and is enclosed in a blue rounded rectangular border.

SmitCare

PROMOTING INDEPENDENCE
FOR DAILY LIVING

WWW.SMITCARE.CO.UK




SmitCare

Welcome

to the SmitCare Trade Catalogue

incorporating





In the following pages you will find a range of products designed to assist people in their daily lives and to offer solutions for carers and professional advisors.

If you are unsure of the suitability of equipment, please contact us for advice.



CONTENTS

Bedroom	2-5
Furniture Raisers	6-10
Bathroom & Toileting	11-16
Mobility	17-19
Grab Rails	20
Moving & Handling	21-25

Bedroom

DreamMaster Light with Optional Knee-Break and Side Rails



DreamMaster Mattress Elevator fitted to a divan bed, with a Left Side Rail

The award winning DreamMaster raises the standard for mattress elevators and allows greater independence for people in their own home.



The unique features of the DreamMaster are retained:

- Fits single and double beds up to 150cm (5') wide
- Can be fitted to a wide range of domestic beds, including divan and slatted bed bases
- Optional side rails raise with the backrest
- Optional Knee-Break to prevent user moving down the bed
- Market leading elevation range, from 5° to 85°, to allow it to be almost flat.
- Virtually silent in use
- Maximum weight rating of 40 stone (254kg)

Code	Description	Price Each
DM300	DreamMaster Light	£

DreamMaster Light Knee-Break

The movement of the Knee-Break is controlled independently from the DreamMaster using the four-button hand control supplied.

The Knee-Break has a maximum vertical height of 20cm. This product can only be used in conjunction with the DM300 Mattress Elevator.



Code	Description	Price Each
DM400	DreamMaster Knee-Break	£

DreamMaster Light Side Rails




The Left and Right Side Rails can be fitted to the DreamMaster Light Mattress Elevator and raise with the backrest to provide a firm support for repositioning and sit to stand transfers. The side is identified from the user's position lying on the bed.

Code	Description	Price Each
DM300/L	DreamMaster Light Left Rail	£
DM300/R	DreamMaster Light Right Rail	£

DreamMaster Light Mattress Support Bars

Offer more support to softer mattresses



Code	Description	Price Each
DM300/MSB	DreamMaster Light Mattress support bars	£

DreamMaster Light spares

Code	Description	Price Each
DM300/HC2	Hand Control for DM300	£
DM300/TB	Transformer for DM300	£
DM300/ACT	Actuator for DM300	£
DM300/HC4	Four Button Hand Control for DM400	£
DM400/ACT	Actuator for DM400	£
DM300/BS	2.5m brown straps (pack of ten)	£
DM300/CS	3.5m cream straps (pack of ten)	£
DM300/CTS	Cream Tail Straps (pack of ten pairs)	£
DM300/BKT	Pair of Strap Brackets with straps attached	£
DM300/STRAP	Complete set of straps & brackets for DM300	£

Bedroom

SmitCare Backrest

Secure and comfortable backrest which can be adjusted from a low lying angle to an upright sitting position. Ideal for Health and Social Services.

The back rest can be set at 30, 40, 50, 60 or 70 degrees.



Code	Description	Price Each
SAC830	SmitCare Backrest	£

SmitCare Bed Bar



A practical supporting grab bar to benefit the frail and unstable. The bar height is adjustable for the comfort and stability of the user.

The large wooden base slides between the bed base and mattress. The inverted 'T' safety bar eliminates the risk of entrapment. Supplied with fitting straps which locate around the bed base securing the Bed Bar to the bed.

The Bed Bar can be set at either 36cm or 45cm above the height of the bed base.



Code	Description	Price Each
SBH1/ST	Bed Bar with straps	£

Medeci Pillow Lifter Light

The Medeci Pillow Lifter Light helps the user to sit up in bed. It is quiet, comfortable and stable in any position on the up/down cycle.



The mattress pad has a vapour permeable, four-way stretch waterproof cover, which secures to the frame and a fitted polycotton sheet with a sewn-in pillow pocket.

The aluminium frame, mattress and covers weigh 13kg. The strap may be used to attach the PillowLifter Light to the mattress/bed base if additional security is required. The strap is long enough to secure the PillowLifter to a 150cm bed. An additional strap can be used to fit to larger beds.

Backrest size: 750mm wide, 720mm tall
Apron size: 750mm wide, 280mm long
Angle of operation: 5° to 65°

Spares are available on request.

**WITHOUT
FITTED SHEET
IN PLACE**



Code	Description	Price Each
MPL100	Medeci PillowLifter Light	£

Medeci Pillow Lifter Light spares

Code	Description	Price Each
MPL100/FS	Fitted Sheet	£
MPL100/WC	Waterproof Cover	£
MPL100/ST	Spare Strap	£
MPL100/HC	Hand Control	£
MPL100/TB	Transformer Box	£
MPL100/ACT	Actuator	£

Furniture Raisers

Medeci Furniture Raisers

- Suitable for a wide range of castors and legs
- Adaptors included
- Integral "twist" height adjustment over 3" range
- Reference guide on inner barrel to ensure easy height adjustment
- Base flanges designed to allow furniture to be moved back into position after fitting of raisers



**EASY
HEIGHT
ADJUSTMENT
GAUGE**




The Medeci MFR1 raisers allow a raise from 2" (50mm) to 5" (127mm) for furniture supported on the base of the raiser cup, and 4" (100mm) to 7" (177mm) less the height of the existing castor for divans and armchairs; most castors are about 1.75" (45mm) in height. The raisers are suitable for all traditional and modern cleaning systems.

Code	Description	Price Each
MFR1	Medeci Furniture Raisers (set of 4) inc. adaptors	£
MUCA	Medeci Universal Castor Adaptors (set of 4)	£

Medeci Link Arms

The MLA1/U Universal Link Arm assembly for use with the Medeci Furniture Raisers and Clip & Stack Raisers.

The Link Arm assembly has been developed to enhance the versatility of the popular Medeci Furniture Raisers. The aluminium arms are strong, light and can be machine cleaned. The raisers simply clip in, without the need for tools.



The Link Arm assembly imparts additional stability and security. The MLA1/U can be used in all situations where the MFR1 or CFR2 are suitable and there is no requirement for clear access under the furniture.

The standard unit is supplied in a chair formation as illustrated but additional short and long stretcher bars are available, with nuts and bolts, to allow use in a wide range of formations and applications.

**MEDECI
FURNITURE
RAISERS
SOLD
SEPARATELY**



Code	Description	Price Each
MLA1/U	Medeci Link Arm Universal	£
MLA1/S	40cm Stretcher Bar (pack of 6)	£
MLA1/L	100cm Stretcher Bar (pack of 6)	£
MLA1/WS	Wingnut & bolt set for MLA1/U	£

Furniture Raisers


Medeci Clip & Stack

The Medeci Clip & Stack furniture raisers are designed to fit furniture with large diameter legs or with flat slider plates. They attach directly to the underside of the leg or floor plate, providing a neat, cost effective and unobtrusive raise for furniture that does not fit other raisers.

The top disc attaches to the furniture using the strong adhesive pad, or it can also accept a screw fixing if desired. Once the top disc is in place, the clip-on sections are added to create a stack of the required height, with a slip resistant boot supplied at the bottom for use on hard flooring.

The Clip & Stack Raisers can also be linked with the Medeci Link Arms for additional stability and security.

Height options - 40mm, 65mm or 90mm)




Code	Description	Price Each
CFR2	Clip & Stack Raisers (set of 4)	£
CFR2/TS	Spare Top Disc (set of 4)	£

Multi-purpose Raisers

Raises beds, chairs and settees

The Langham Multi-Purpose Raisers are the modern descendants of the original Langham range, designed to fit chairs, settees and beds. Used as a modular system building from four separate elements, the Langham Multi-Purpose Raisers replace eleven original dedicated Langham products.



- (A) MPR Unit
- (B) Corner Support Well
- (C) Pivoting Corner Support
- (D) Clip-on Height Sections
- (E) Spreader Link Bar




Code	Description	Price Each
LPMPR1	Set of 2 Linked Supports	£
LPMPR2	6 x 20" Spreader Bars	£
LPMPR3	12 x Clip-on Height Adjusters	£
LPMPR4	6 x 42" Spreader Bars	£

Elephant Feet

Gordon Ellis Leg Raising Cones. Grey in colour to compliment most furnishings and made from strong recyclable polypropylene. Centres chair, table, bed or settee legs in a safe position, though with this system we recommend beds should be against at least two walls.


Two heights offered:
90mm and 140mm (3 1/2" and 5 1/2"). Sold in sets of 4.



Code	Description	Price Each
SE6070	Elephant Feet 90mm (3 1/2") (carton of 4)	£
SE6074	Elephant Feet 140mm (5 1/2") (carton of 4)	£

Langham Chair Raiser

Langham's original chair raiser for chairs with legs. Supplied with height adjustment blocks.



Code	Description	Price Each
LLR1	Langham Chair Raiser Adjustable Height	£

Furniture Raisers

Langham Bed Raiser

Langham's original bed raiser for beds with legs. Supplied with height adjustment blocks. Two pairs of linked supports per box.



Code	Description	Price Each
LLR3	Bed Raising Unit Adjustable Height (for bed with legs)	£

Langham Grip-On Raisers

Simple to fit individual chair and table leg raisers. They can be used time and time again and once fitted hold firmly to the leg until deliberately removed. When fitted correctly they will cause no damage to the furniture. Height is adjusted simply by adding or removing sections of the raiser to achieve a 35, 60 or 85mm raise. This range is tested safe to 318kgs (50 stones).



Code	Description	Price Each
LLR100	38mm (1 1/2") Leg. dia.	£
LLR200	50mm (2") Leg. dia.	£
LLR300	64mm (2 1/2") Leg. dia.	£

Myco Ultra Swivel Bather®



The Myco Swivel Bather® range has been enhanced by the development of two products. The new Ultra Swivel Bather is a result of the popularity of the existing adjustable models. The Ultra Swivel Bather is width adjustable between 26" and 30", has a weight rating of 160kg (25 stone) but differs from the existing models because the frame sections are made from aluminium for increased ease of recycling and durability. The seat has also been restyled for comfort while retaining drainage holes and complying with regulations.

The Ultra Swivel Bather has two parts; a base frame that rests on the bath and a seat unit that rotates on the base. The two parts separate easily for handling and storage. The total weight of the product is 5kg.

The Ultra Swivel Bather is supplied with the locking lever located on the right hand side of the seat unit and locks in both the transfer and bathing positions. The locking lever can easily be changed to the left side by removing a bolt, changing the lever position and replacing the bolt.

The Swivel Bather® fits on any standard bath, no fixtures or fittings are required and the unit is easily removed when not required.

Dimensions:

- Seat width: 445mm (17½")
- Backrest height: 355mm (14")
- Seat depth: 405mm (16")
- Width between arms: 535mm (21")
- Seat height above side of bath: 55mm (2.2")
- Distance between arms of base frame (measured at outside edge of bath): 535mm (21")
- Max user weight rating: 160kg (25st)



Code	Description	Price Each
MUSB	Myco Ultra Swivel Bather	£

SmitCare Ultra Bath Step

The Ultra Bath Step is comprised of four 25mm (1") sections that push together to form a step up to 100mm (4") high. Further sections can be added up to a maximum of 200mm (8"). There are indents on each section to make it easier to separate the sections.

The upper surface of the sections is textured to provide a firm footing and the bottom section has rubber feet for slip-resistance.

Dimensions: 350 x 450mm (14 x 18")



Code	Description	Price Each
SUBS	SmitCare Ultra Bath Step	£

Myco Adjustable Bathboard



Myco has developed this high quality, elegantly simple adjustable bathboard, for use on baths from 660 to 760mm (26" to 30") wide.

Simply adjust the length of the board to suit the width of the bath, insert the handle to secure, and then adjust the two brackets to fit the internal walls of the bath.

The bathboard is strong, comfortable and is suitable to fit to most domestic baths.

The bathboard is 300mm wide.

Note: The handle is included with each board.



Code	Description	Price Each
MBAX	Adjustable Bathboard	£
MBAX/HAN	Myco Bathboard Handle	£

Bathroom & Toileting

Yateley Shower Stool

The Yateley Shower Stool is designed for use in a shower cubicle. The stool features drainage holes in the seat, height adjustable legs and slip-resistant rubber ferrules. The seat is removable for cleaning.

Dimensions:

Seat size: 37x28cm
 Width between arms: 47cm (18½")
 Height adjustment: 50-64cm (19-25")
 Max footprint 58 x 42cm (23" x 16½").



Code	Description	Price Each
SW127/A	Yateley Shower Stool	£

Medeci Ultra Bathseat

The Medeci Bathseat has always been popular and successful because of the width adjustable seating area. Although the suckers give a secure fitting to the bath, the ability to expand the seat to fit the internal width of the bath provides the user with an enhanced sense of security and stability.

Previously, to adjust the width of the seat, two adjusters had to be rotated simultaneously. The revised product has been simplified and has just one adjuster, making fitting easier and quicker.

The Medeci Ultra Bathseat is available in three heights to meet the most popular demands. The Bathseat is supplied part assembled to allow a reduction in packaging materials and storage space, and simply requires the legs to be inserted before use.



Adjustable width:
 480 - 640mm (19" to 25")
Seat depth:
 220mm (8½")



Code	Description	Price Each
MUBS6	6" Medeci Ultra Bathseat (150mm)	£
MUBS8	8" Medeci Ultra Bathseat (200mm)	£
MUBS12	12" Medeci Ultra Bathseat (355mm)	£

Bathroom & Toileting

Aquarian Bathboard & Handle

The Aquarian Bathboard is made from strong extruded plastic tube and strengthened with internal plywood slats. This comfortable easy to clean bathboard is a most popular choice. Deep, easy to fit cleats.

Dimensions:

Width: 230mm (9")

Available in three Lengths: 660mm (26"), 685mm (27"), 710mm (28").



Code	Description	Price Each
ASB6226	Slatted Bathboard 660mm (26")	£
ASB6227	Slatted Bathboard 685mm (27")	£
ASB6228	Slatted Bathboard 710mm (28")	£

Aquarian Showerboard

Produced in the same strong hygienic materials as the Aquarian Bathboard, the Aquarian Showerboard fits comfortably into most baths. It is easy to clean and gives an ideal seating platform for those requiring a larger seat area.

Dimensions:

Width: 350mm (14")

Available in three lengths: 660mm (26"), 685mm (27"), 710mm (28")



Code	Description	Price Each
AQSB6626	6-Slat Showerboard 26"	£
AQSB6627	6-Slat Showerboard 27"	£
AQSB6628	6-Slat Showerboard 28"	£

Ashby Quick-Fit Contract Raised Toilet Seat

The Ashby Quick-Fit toilet seat has been developed with easy and effective cleaning in mind. Designed to make fitting even easier; the tightening handwheels are smaller than those fitted to the Ashby Super Contract. With the integral front lip it offers a secure fit to almost all shapes of toilet pan. Produced in re-cyclable polyethylene and impregnated with anti-bacterial properties. All round, this RTS represents great value to the equipment provider, especially those re-issuing cleaned products. Available in three sizes, 50mm (2"), 100mm (4") and 150mm (6").



Code	Description	Price Each
AQFC2	50mm (2")	£
AQFC4	100mm (4")	£
AQFC6	150mm (6")	£

Ashby Contract Raised Toilet Seat

- Comfortable to use and warm to touch
- Produced in easy clean Polyethylene
- Impregnated with Anti-Bacterial Irgasan® (Kills MRSA)
- Available in 50mm, 100mm and 150mm (2", 4" and 6")
- Safe reliable fixings
- Sealed unit, will not take in water during cleaning
- Autoclavable for one minute at 90°C
- Deep frontal access for personal hygiene
- Fits most standard toilet bowls
- Tested safe to 190kgs (30 stones)
- Three year manufacturers guarantee
- Contract version in opaque white



Super Ashby Contract: Easy Screw Fix

Code	Description	Price Each
ATSC2E	50mm (2")	£
ATSC4E	100mm (4")	£
ATSC6E	150mm (6")	£

Super Ashby Contract: Push on Fix (not illustrated)

Code	Description	Price Each
ATSC2N	50mm (2")	£
ATSC4N	100mm (4")	£
ATSC6N	150mm (6")	£

This extremely popular market leader is now produced in Irgasan® impregnated Anti-Bacterial polyethylene and the Easy-Screw clip covers are moulded on to prevent them being lost when the seat is removed. The Ashby Contract is specifically aimed at the Hospital / Social Services market and offered in an opaque pearl finish.

Bathroom & Toileting

Apex Mobile Commode

The Apex Mobile Commode is a strong, reliable product for use as a static commode, or to transfer somebody to the bathroom for over-toilet use. The armrests are removable to allow sideways transfer, the footrests are conventional swing-away & lift-off, the teardrop shaped aperture allows better access and the large round potty has a secure lid.



Dimensions:

Seat height: 500mm (19.5")

Seat size: 435mm wide x 460mm deep (17" x 18")

Footprint: 460mm x 510mm (18"x20")


Aperture: 200mm diameter, opening to 280mm front to back.

Four brake castors, diameter 130mm (5")



Code	Description	Price Each
SAC360	Apex Mobile Commode	£

Metal Half Step



Anti-slip tread plate, non-rust, In two standard sizes. Height 115mm (4 1/2").



Code	Description	Price Each
LLR9	Standard Step 405 x 760mm (16" x 30")	£
LLR10	Large Step 610 x 760mm (24" x 30")	£

Derby Half Step




Polypropylene with strong, non-slip rubber mat. Adjustable height to enable use on uneven surfaces. Height 115mm (4 1/2").



Code	Description	Price Each
SE66091	Derby Half Step 740 x 440mm (29" x 17")	£

Smitcare Threshold Ramps

SmitCare Threshold Ramps are manufactured in anodised aluminium and trimmed at the edges with a coarse finish strip to the side up stands preventing the risk of slippage. The treads are slotted to facilitate drainage and stop the build up of ice. Easily transportable in the back of a car. Threshold ramps are available in the following sizes.



Dimensions:	Length	Width	Weight
	2ft (610mm)	30" (76cm)	4.8kg
	3ft (915mm)	30" (76cm)	7.1kg

Code	Description	Price Each
TR2	2ft Threshold Ramp	£
TR3	3ft Threshold Ramp	£

Mobility

Smitcare Folding Ramps

SmitCare Folding Ramps are manufactured in anodised aluminium and trimmed at the edges with a coarse finish strip to the side up stands preventing the risk of slippage. The treads are slotted to facilitate drainage and stop the build up of ice. Supplied with handle and lock.

Easily transportable in the back of a car. Folding ramps are available in the sizes opposite.

Ramp width 29" (740mm) when open.



Dimensions:	Size	Weight
	3ft(915mm)	9.2kgs / 20lbs
	4ft(1220mm)	12.2kgs / 27lbs
	5ft(1525mm)	15.2kgs / 34lbs

Code	Description	Price Each
FWR3	3ft Folding Ramp	£
FWR4	4ft Folding Ramp	£
FWR5	5ft Folding Ramp	£

Ramp Wedge



Portable freestanding ramp for use at low door thresholds.

Code	Description	Price each
LLR12/10	1" (25mm) Ramp Wedge 30"	£
LLR12/20	2" (51mm) Ramp Wedge 30"	£
LLR12/30	3" (76mm) Ramp Wedge 30"	£

Tri-Wheel Partner

This foldable, lightweight walker can be supplied with a hook on basket or bag. Height adjustable handles and locking handbrakes. Large microcellular tyres which do not need to be pumped up. Folds flat to fit into any car boot.

Height adjustment of handle 825mm - 914mm (32½"-36"). Weight 6.4kgs (14lbs).



Code	Description	Price Each
SAC4	Tri-Wheel Partner	£
SAC2/BAS	Basket for SAC2	£
SAC6	Tri-Wheel Partner with Bag	£

Lightweight Rollator

The lightweight rollator benefits from an aluminium frame and 8" (20cm) castors, ergonomic handles, lock-on lever brakes, padded removable back rest, padded seat, and an under seat wallet and bag.



Dimensions:

Height adjustable handles: 32½" - 37" (82.5cm - 94cm).
 Overall width: 24" (61cm).
 Width between arms: 20½" (52cm).
 Seat height: 21" (52.5cm).
 Weight: 7kgs (15.4lbs).



Code	Description	Price Each
LSR24	Lightweight Rollator	£

Buckingham Caddy

The Buckingham Caddy was designed as an aid to safely carry a dinner plate, or cup using the supplied attachment. There is an ample storage area underneath for personal items. Fits most types of rigid frame.



Centre of gravity remains behind the front wheels to maintain user stability



Code	Description	Price Each
BWFC1/U	Walking Frame Caddy	£

Grab Rails & Supports

Aquarian Fluted Grab Rails with Shroud

A very strong and comfortable fluted / ribbed handrail with a 33mm (1 1/4") diameter handle grip. The 70mm (2 3/4") rose is countersunk with three 6mm holes. This is covered after fitting by a secure press-down screw cover, referred to as a shroud. Available in three lengths. Supplied in boxes of 25 rails per size. Not recommended for exterior use.



Code	Description	Price Each
AQGR12	Shrouded Plastic Fluted Handrail 305mm (12")	£
AQGR18	Shrouded Plastic Fluted Handrail 460mm (18")	£
AQGR24	Shrouded Plastic Fluted Handrail 610mm (24")	£


Flat End Grab Rails

Coated steel grab rails with flat ends with two screw holes to secure to the wall.

There is often some discrepancy in the description of these rails, we quote the overall length of the rails, including the flat end sections. For example, the MGR1A – 18" Metal Grab Rail is the overall length of the rail. The approximate grab length is about 12" for this rail.



		Price Each
MGR1/A	Metal Grab Rail 460mm (18")	£
MGR2/A	Metal Grab Rail 610mm (24")	
MGR3/A	Metal Grab Rail 700mm (28")	



Locomotor Ultra Slide Sheets

The ULTRA range of slide sheets is available in seven sizes and allows three basic patient handling tasks without the need to lift and without the risk of sheering the patient's skin. The multidirectional glide and ultra low friction of the sheets makes sliding transfers extremely easy to perform. Can be used with sliding boards to ease a sliding transfer from chair to wheelchair etc or to aid lateral transfer from bed to trolley using a bridging board.

All Sliding Sheets can be autoclaved or washed above 160°F.

Benefits:

- Turning a patient: for X-Ray, surgical procedures, recovery, pressure care, change of dressing or cleaning.
- Sitting a patient up on a bed: for feeding, relief of chest pain etc...
- Moving a patient on or off a bed, trolley or examination couch



Code	Description	Price Each
LOCO-001X	Standard 122cm x 71cm	£
LOCO-002X	Long 145cm x 71cm	£
LOCO-003X	Wide 122cm x 100cm	£
LOCO-012	Compact 72cm x 70cm	£
LOCO-013	Maxi 195cm x 70cm	£
LOCO-067	Maxi Plus 195cm x 140cm	£
LOCO-076	Mini 60cm x 41cm	£

Locomotor Quilted Uni-Slide

Available in two sizes, the Uni Slide prevents patients from sliding down in a chair or bed. The tubular Uni Slide allows an assisted slide back in to the chair or back up the bed, and then locks to prevent the patient sliding forward.



Code	Description	Price Each
LOCO-190	(45cm x 40cm) Quilted Uni Slide Chair	£
LOCO-191	(80cm x 70cm) Quilted Uni Slide Bed	£

Locomotor Quilted Flat Uni Slides

Developed as an alternative to the tubular Uni Slides, the quilted Flat Uni Slide is available in two sizes and incorporates a series of handholds which allow a positive grip when assisting the patient.



Code	Description	Price Each
LOCO-192	(90cm x 40cm) Quilted Flat Uni Slide Chair	£
LOCO-193	(160cm x 70cm) Quilted Flat Uni Slide Bed	£

All Select Healthcare Moving & Handling products are available from SmitCare

Moving & Handling

Rotastand

The Rota Stand has been developed to allow an assisted transfer from one seated position to another. The frame shin pad and turn-table combine to give optimum safety to both carer and patient during the transfer. The central wheels are a unique innovation as they allow the Rota Stand to be rolled up to the user even if their feet are already in position.

Features: Height adjustable, foot pedal to secure turn-table, central wheels for easy positioning and adjustable shin-pad for correct stand position.



Code	Description	Price Each
LOCO-102	Rotastand	£

Rotastand Compact



The new Rotastand Compact is a lightweight, folding version of the Rotastand range. Weighing just 10kg and with a handle that folds flat, it is much easier to transport and store.

As with previous models, the shin pad and handle are height adjustable to suit the user, the wheels are now covered and the stabilizing bar is retained to allow a one person transfer.



Code	Description	Price Each
LOCO-600	Rotastand Compact	£

Locomotor Handling Belts & Slings



Locomotor have developed this extensive range of handling slings and belts to cover most of the more commonplace requirements in patient handling. Should you require more details on this range we can supply a dedicated Locomotor brochure by return.

Code	Description	Price Each
LOCO-006	(A) Maxi Belt 24" - 45" (600-1145mm)	£
LOCO-077	(B) Maxi Plus Belt 36" - 60" (915-1525mm)	£
LOCO-098	(C) Maxi Deluxe Belt 24" - 45" (610-1145mm)	£
LOCO-007	(D) Mini Belt 20" - 45" (508-1145mm)	£
LOCO-114	(E) Mini Deluxe Belt 20" - 45" (508-1145mm)	£
LOCO-008	(F) Handling Sling 62 x 22cm	£
LOCO-033	(G) Economy Belt up to 45" (1145mm)	£

Boomerang Transfer Boards

Top quality plywood transfer boards, safe to 190 kilos (30 stones) when used correctly. Feathered at each end and trimmed with a non-slip rubber strip. Boards have a low friction surface. Economically priced, packed in boxes of 10.

Dimensions:

Length overall 740mm (29"), Width 260mm (10").




Code	Description	Price Each
STBTB	Boomerang Transfer Board Boxed in 10s	£

Moving & Handling

Straight Transfer Boards

The Transfer Boards are straight boards, manufactured from top quality plywood. They are available in two lengths for suitability and have feathered edges for ease of use.

The longer 76cm board also features cut out hand-holds.



Code	Description	Price Each
DTB24S	Transfer Board 610 x 203mm [24" x 8"]	£
DTB30S	Transfer Board 750 x 203mm [29 1/2" x 8"]	£

Locomotor Duo Transfer Slide


For angled or straight transfer from one seated position to another, very smooth and virtually unbreakable. Low friction surface for easy smooth slide.



Code	Description	Price Each
LOCO-401B	Loco Duo Transfer Slide	£
LOCO-400B	Duo Transfer Slide with Handles	£

Turntables


Available in two sizes (with a plastic or fabric surface), the turntable allows clients to be pivoted safely without the need to move their feet. The turntable rotates smoothly and safely with minimal effort, which reduces the risk of loss of balance and injury.



Code	Description	Price Each
LOCO-034	Fabric turntable - 30cm	£
LOCO-035	Fabric turntable - 38cm	£
LOCO-038	Plastic Turntable - 30cm	£
LOCO-039	Plastic Turntable - 38cm	£

Rota Cushions

A soft, padded, flexible turning disc for use on car seats, chairs and beds. The Rota Cushion allows a smooth swivel transfer without painful jarring or twisting. The central swivel joint attaches to the sliding discs with velcro, so that they can be removed for washing. Available in two sizes.



Code	Description	Price Each
LOCO-068	Rota Cushion - 38cm	£
LOCO-078	Rota Cushion - 45cm	£

SUPA Slide Sheets

The SUPA Slide sheets have been developed to help with the increasing necessity for infection control. The SUPA Slide has the same performance as the Ultra Slide sheets and permits multiple use with one patient but can be disposed of when it becomes soiled. The sheets are supplied in a dispenser box for easy removal and can either be used in pairs with one on top of the other, or folded in half to form a tube. Two sizes are available.



Code	Description	Price Each
LOCO-200	Box of 100 SUPA Slide - 200x75cm	£
LOCO-201	Box of 100 SUPA Slide Wide - 200x100cm	£

Terms and Conditions of Trading

Minimum Order Value: £10 (nett goods value)

Pricing and Specifications:

All Goods in this catalogue are subject to being unsold and alteration in price or specification without prior notice. Prices in this catalogue refer to single units. In some instances the units are only supplied at this price in pre-packaged dispatch quantities.

Value Added Tax is not included in the price list and will be charged at the appropriate rate on invoice (goods plus carriage and handling). Certain items are allowable for zero VAT rating. Charities may benefit from this provided the relevant VAT Zero Rate Certificate is completed and the goods are allowable for zero rating by H.M. Customs and Excise. Certificates are available from SmitCare on request.

Measurements and Weights:

Where stated are approximate and are not meant to be exact.

Copyrights where applicable are shown.

Carriage & handling: Orders of £150 or over (nett value) carriage paid U.K. mainland only except postcodes listed in Appendix 1 below, where carriage charges attract 100% surcharge. Orders to U.K. mainland with a goods value below £150 will be subject to an £8.50 charge (Appendix 1 surcharge applies) to cover administration, packing and carriage. All deliveries to Northern Ireland will be subject to an £8 consignment charge; orders of £150 or over will be carriage paid, with orders under £150 charged at cost, with a minimum of £8.50. It is an H.M Customs requirement that VAT is charged on Carriage & Handling. All other countries and off-shore islands, including Channel Islands and Isle of Man, ex-warehouse at cost.

Forwarding: By a National Carrier or by Post at our discretion. Next day service available at £14.50 for 1st carton plus £4.50 per carton thereafter U.K.

Mainland only (see Appendix 1 opposite). N.B. This charge is applicable irrespective of order value and should be arranged by 12 noon on day of despatch. Goods are normally despatched within 5 working days from receipt of order and delivered within the next 5 days.

Payment: Approved credit accounts 30 days from date of invoice.

Appendix 1: Postcodes where 100% surcharge applies: HS1-9, IM, IV1-49, IV51-56, KA27-28, KW1-3, KW5-17, PA20-49, PA60-78, PH17-26, PH30-44, TR21-25, ZE1. Carriage paid if order value over £300.

Receipt of Goods:

All packages must be checked on arrival for damage and appropriately signed for. Should there be cause to raise a claim our offices must be notified within 24 working hours of the delivery. Carriers will not consider claims outside these parameters. Should examination not be possible on receipt, goods should be signed for unexamined or refused.

All Claims must be notified by telephone and confirmed in writing within 3 days. This is vitally important as any claim is made on the purchasers' behalf.

Booked Deliveries If customers require deliveries booked-in there will be a charge levied in addition to applicable carriage charges at carriers cost.

Damage: Any parcels must be signed for as above. If you are unable to check the contents on delivery the consignment must be signed for as 'unexamined' and any damage notified within 24 hours.

Discrepancy: Shortages on delivery must be notified within 3 days of delivery.

Loss: Non-delivery of goods must be notified within 10 days of Invoice /Advice Note date.

Returns: Please telephone our Customer Services department if you wish to return goods. We reserve the right to charge a handling fee of 10% on goods ordered in error. Goods will only be eligible for credit if they are in boxed, resaleable condition. Return carriage costs of incorrectly ordered goods is at the expense of the customer.

SmitCare is an Equal Opportunities employer.

Hours of Business: Monday to Friday 9am - 5pm





Products to Assist Quality Care in our Community

SmitCare Ltd, 14 Camberley Business Centre, Camberley, Surrey, GU15 3DP

Telephone: 01276 677888

email: sales@smitcare.co.uk

WWW.SMITCARE.CO.UK